





DEWIMED®



CONTAINER SYSTEM


DEWISTERI

DEWISTERI CONTAINER	3-10
DEWISTERI CONTAINER	
DEWISTERI SISTEMAS DE CONTENEDORES	

DEWIMED CONTAINER	11-19
DEWIMEDCONTAINER	
DEWIMED CONTENEDORES	

SIEBKÖRBE	19-28
WIRE BASKETS AND TRAYS	
BANDEJAS DE ACERO INOXIDABLE	

ZUBEHÖR.....	29-32
ACCESORIES	
ACCESORIOS	



DEWISTERI

Safe, Saver, Savest

Seguro, Más Seguro, Lo Más Seguro

Sicher, Sicherer, am Sichersten

Our standard model offers a simple set-up, easy handling and an optimal price/performance ratio.

Based on these values a high end model was developed. With its safety lid a maximum of protection during transport and storage can be guaranteed. Furthermore it reduces the chance of contamination significantly.

El modelo estandar ofrece una estructura simple, el máximo de facilidad en el uso y un óptimo de la relación precio-calidad.

En base de estos valores un modelo de tecnología alta fue desarrollado. Con su tapa de seguridad puede garantizar un máximo de protección para el transporte así como el almacenamiento. Además reduce el riesgo de contaminación significativamente.

Das Standardmodell kann durch den einfachen Aufbau, maximale Bedienerfreundlichkeit und ein optimales PreisLeistungsverhältnis überzeugen. Darauf aufbauend wurde ein High-End-Modell entwickelt, dass durch den Sicherheitsdeckel nicht nur zusätzlichen Schutz bei Transport und Lagerung bietet, sondern auch die Kontaminierungswahrscheinlichkeit deutlich reduziert.

Standard Model
Modelo Estándar
Standardmodell



Safety model
Modelo de seguridad
Sicherheitsmodell




DEWISTERI




Security Handle Mango de Seguridad Sicherheitshandgriff

- Extra long handle with silicone tube
- Protects against burning and provides more grip
- Mango más largo con silicona
- Protege de quemadura y preventa el deslizamiento
- Extra langer Griff mit Silikonüberzug
- Schützt vor Abrutschen und Verbrennungen




Intuitive Handling Uso intuitivo Intuitive Bedienung

- Self-explaining change of filter
- Color codes for easy administration
- Fácil cambio del filtro
- Código de colores por una administración transparente
- Selbsterklärender Filterwechsel
- Farbcodes für einfache Verwaltung



Aluminium Aluminio Aluminium

- Optimal thermal conductivity and drying process
- Proven durability
- Conductividad termica y proceso de secamiento óptimál
- Durabilidad probado
- Optimale Wärmeleitfähigkeit und Nachtrocknung
- Bewährte Beständigkeit




Hygiene hygiene Hiegiene

- Easy to clean
- Smooth surfaces
- Completely detachable
- Fácil para limpiar
- Superficios planos
- Desmontable completamente
- Leicht zu reinigen
- Glatte Flächen
- Komplett zerlegbar

DEWISTERI

Aluminium 1/1 Container Aluminium 1/1 Container Contenedor 1/1 de Aluminio

Boden ungelocht
Bottom non perforated
Fondo sin perforar



5200 **Zyklen**/Cycles/Ciclos

100	135	150	200	260
außen / outside / exterior 580 x 280 x 100 mm	außen / outside / exterior 580 x 280 x 135 mm	außen / outside / exterior 580 x 280 x 150 mm	außen / outside / exterior 580 x 280 x 200 mm	außen / outside / exterior 580 x 280 x 260 mm
innen / inside / interior 550 x 265 x 84 mm	innen / inside / interior 550 x 265 x 115 mm	innen / inside / interior 550 x 265 x 130 mm	innen / inside / interior 550 x 265 x 180 mm	innen / inside / interior 550 x 265 x 240 mm

PTFE System / PTFE System / Sistema de PTFE

	88.510.10DAP		88.510.13DAP		88.510.15DAP		88.510.20DAP		88.510.26DAP
	88.520.10DAP		88.520.13DAP		88.520.15DAP		88.520.20DAP		88.520.26DAP
	88.530.10DAP		88.530.13DAP		88.530.15DAP		88.530.20DAP		88.530.26DAP
	88.540.10DAP		88.540.13DAP		88.540.15DAP		88.540.20DAP		88.540.26DAP
	88.550.10DAP		88.550.13DAP		88.550.15DAP		88.550.20DAP		88.550.26DAP
	88.560.10DAP		88.560.13DAP		88.560.15DAP		88.560.20DAP		88.560.26DAP

Barriere System / Barrier System / Sistema de Barrera

	88.510.10DAB		88.510.13DAB		88.510.15DAB		88.510.20DAB		88.510.26DAB
	88.520.10DAB		88.520.13DAB		88.520.15DAB		88.520.20DAB		88.520.26DAB
	88.530.10DAB		88.530.13DAB		88.530.15DAB		88.530.20DAB		88.530.26DAB
	88.540.10DAB		88.540.13DAB		88.540.15DAB		88.540.20DAB		88.540.26DAB
	88.550.10DAB		88.550.13DAB		88.550.15DAB		88.550.20DAB		88.550.26DAB
	88.560.10DAB		88.560.13DAB		88.560.15DAB		88.560.20DAB		88.560.26DAB

Ventil System / Valve System / Sistema de Válvula

	88.510.10DAV		88.510.13DAV		88.510.15DAV		88.510.20DAV		88.510.26DAV
	88.520.10DAV		88.520.13DAV		88.520.15DAV		88.520.20DAV		88.520.26DAV
	88.530.10DAV		88.530.13DAV		88.530.15DAV		88.530.20DAV		88.530.26DAV
	88.540.10DAV		88.540.13DAV		88.540.15DAV		88.540.20DAV		88.540.26DAV
	88.550.10DAV		88.550.13DAV		88.550.15DAV		88.550.20DAV		88.550.26DAV
	88.560.10DAV		88.560.13DAV		88.560.15DAV		88.560.20DAV		88.560.26DAV

Zubehör
Accessories
Accesorios

PTFE Filter/PTFE Filter/Filtro de PTFE

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenancance free/Sin Mantenimiento

Art. No. 88.701.02



Barriere Filter/Barrier Filter/Filtro de Barrera

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenancance free/Sin Mantenimiento

Art. No. 88.701.03




Ventil Filter/Valve Filter/Filtro de Válvula

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenancance free/Sin Mantenimiento

Art. No. 88.701.04



Container Deckel
Container Lid
Tapa de Contenedor



PTFE System PTFE System Sistema de PTFE		Barriere System Barrier System Sistema de Barrera		Ventil System Valve System Sistema de Válvula	
	88.564.01DAP		88.564.01DAB		88.564.01DAV
	88.564.02DAP		88.564.02DAB		88.564.02DAV
	88.564.03DAP		88.564.03DAB		88.564.03DAV
	88.564.04DAP		88.564.04DAB		88.564.04DAV
	88.564.05DAP		88.564.05DAB		88.564.05DAV
	88.564.06DAP		88.564.06DAB		88.564.06DAV

DEWISTERI

Aluminium 3/4 Container Aluminium 3/4 Container Contenedor 3/4 de Aluminio

Boden ungelocht
Bottom non perforated
Fondo sin perforar



5200 **Zyklen**/Cycles/Ciclos

100	135	150
außen / outside / exterior 465 x 280 x 100 mm	außen / outside / exterior 465 x 280 x 135 mm	außen / outside / exterior 465 x 280 x 150 mm
innen / inside / interior 425 x 260 x 84 mm	innen / inside / interior 425 x 260 x 115 mm	innen / inside / interior 425 x 260 x 130 mm

PTFE System / PTFE System / Sistema de PTFE

	88.410.10DAP		88.410.13DAP		88.410.15DAP
	88.420.10DAP		88.420.13DAP		88.420.15DAP
	88.430.10DAP		88.430.13DAP		88.430.15DAP
	88.440.10DAP		88.440.13DAP		88.440.15DAP
	88.450.10DAP		88.450.13DAP		88.450.15DAP
	88.460.10DAP		88.460.13DAP		88.460.15DAP

Barriere System / Barrier Systema / Sistema de Barrera

	88.410.10DAB		88.410.13DAB		88.410.15DAP
	88.420.10DAB		88.420.13DAB		88.420.15DAP
	88.430.10DAB		88.430.13DAB		88.430.15DAP
	88.440.10DAB		88.440.13DAB		88.440.15DAP
	88.450.10DAB		88.450.13DAB		88.450.15DAP
	88.460.10DAB		88.460.13DAB		88.460.15DAP

Ventil System / Valve System / Sistema de Válvula

	88.410.10DAV		88.410.13DAV		88.410.15DAV
	88.420.10DAV		88.420.13DAV		88.420.15DAV
	88.430.10DAV		88.430.13DAV		88.430.15DAV
	88.440.10DAV		88.440.13DAV		88.440.15DAV
	88.450.10DAV		88.450.13DAV		88.450.15DAV
	88.460.10DAV		88.460.13DAV		88.460.15DAV

Zubehör
Accessories
Accesorios

PTFE Filter/PTFE Filter/Filtro de PTFE

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento

Art. No. 88.701.02



Barriere Filter/Barrier Filter/Filtro de Barrera

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento

Art. No. 88.701.03




Ventil Filter/Valve Filter/Filtro de Válvula

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento

Art. No. 88.701.04



Container Deckel
Container Lid
Tapa de Contenedor



PTFE System PTFE System Sistema de PTFE		Barriere System Barier System Sistema de Barrera		Ventil System Valve System Sistema de Válvula	
	88.464.01DAP		88.464.01DAB		88.464.01DAV
	88.464.02DAP		88.464.02DAB		88.464.02DAV
	88.464.03DAP		88.464.03DAB		88.464.03DAV
	88.464.04DAP		88.464.04DAB		88.464.04DAV
	88.464.05DAP		88.464.05DAB		88.464.05DAV
	88.464.06DAP		88.464.06DAB		88.464.06DAV

DEWISTERI

Aluminium 1/2 Container Aluminium 1/2 Container Contenedor 1/2 de Aluminio

Boden ungelocht
Bottom non perforated
Fondo sin perforar



5200 **Zyklen**/Cycles/Ciclos

100	135	150	200	260
außen / outside / exterior 285 x 280 x 100 mm	außen / outside / exterior 285 x 280 x 135 mm	außen / outside / exterior 285 x 280 x 150 mm	außen / outside / exterior 285 x 280 x 200 mm	außen / outside / exterior 285 x 280 x 260 mm
innen / inside / interior 262 x 265 x 82 mm	innen / inside / interior 262 x 265 x 116 mm	innen / inside / interior 262 x 265 x 132 mm	innen / inside / interior 262 x 265 x 182 mm	innen / inside / interior 262 x 265 x 238 mm

PTFE System / PTFE System / Sistema de PTFE

	88.610.10DAP		88.610.13DAP		88.610.15DAP		88.610.20DAP		88.610.26DAP
	88.620.10DAP		88.620.13DAP		88.620.15DAP		88.620.20DAP		88.620.26DAP
	88.630.10DAP		88.630.13DAP		88.630.15DAP		88.630.20DAP		88.630.26DAP
	88.640.10DAP		88.640.13DAP		88.640.15DAP		88.640.20DAP		88.640.26DAP
	88.650.10DAP		88.650.13DAP		88.650.15DAP		88.650.20DAP		88.650.26DAP
	88.660.10DAP		88.660.13DAP		88.660.15DAP		88.660.20DAP		88.660.26DAP

Barriere System / Barrier System / Sistema de Barrera

	88.610.10DAB		88.610.13DAB		88.610.15DAB		88.610.20DAB		88.610.26DAB
	88.620.10DAB		88.620.13DAB		88.620.15DAB		88.620.20DAB		88.620.26DAB
	88.630.10DAB		88.630.13DAB		88.630.15DAB		88.630.20DAB		88.630.26DAB
	88.640.10DAB		88.640.13DAB		88.640.15DAB		88.640.20DAB		88.640.26DAB
	88.650.10DAB		88.650.13DAB		88.650.15DAB		88.650.20DAB		88.650.26DAB
	88.660.10DAB		88.660.13DAB		88.660.15DAB		88.660.20DAB		88.660.26DAB

Ventil System / Valve System / Sistema de Válvula

	88.610.10DAV		88.610.13DAV		88.610.15DAV		88.610.20DAV		88.610.26DAV
	88.620.10DAV		88.620.13DAV		88.620.15DAV		88.620.20DAV		88.620.26DAV
	88.630.10DAV		88.630.13DAV		88.630.15DAV		88.630.20DAV		88.630.26DAV
	88.640.10DAV		88.640.13DAV		88.640.15DAV		88.640.20DAV		88.640.26DAV
	88.650.10DAV		88.650.13DAV		88.650.15DAV		88.650.20DAV		88.650.26DAV
	88.660.10DAV		88.660.13DAV		88.660.15DAV		88.660.20DAV		88.660.26DAV

Zubehör Accessories Accesorios

PTFE Filter/PTFE Filter/Filtro de PTFE

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento

Art. No. 88.701.02



Barriere Filter/Barrier Filter/Filtro de Barrera

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento

Art. No. 88.701.03



Ventil Filter/Valve Filter/Filtro de Válvula

5200 Zyklen/Cycles/Ciclos
Wartungsfrei/Maintenance free/Sin Mantenimiento


Art. No. 88.701.04



Container Deckel
Container Lid
Tapa de Contenedor



PTFE System PTFE System Sistema de PTFE		Barriere System Barrier System Sistema de Barrera		Ventil System Valve System Sistema de Válvula	
	88.655.01DAP		88.655.01DAB		88.655.01DAV
	88.655.02DAP		88.655.02DAB		88.655.02DAV
	88.655.03DAP		88.655.03DAB		88.655.03DAV
	88.655.04DAP		88.655.04DAB		88.655.04DAV
	88.655.05DAP		88.655.05DAB		88.655.05DAV
	88.655.06DAP		88.655.06DAB		88.655.06DAV



DEWIMED Aluminium 1/1 Container

DEWIMED Aluminium 1/1 Container

DEWIMED Contenedor de Aluminio 1/1



1-1200 **Zyklen/Cycles/Ciclos**

100	135	150	200	260
außen / outside / exterior 580 x 280 x 100 mm	außen / outside / exterior 580 x 280 x 135 mm	außen / outside / exterior 580 x 280 x 150 mm	außen / outside / exterior 580 x 280 x 200 mm	außen / outside / exterior 580 x 280 x 260 mm
innen / inside / interior 550 x 265 x 84mm	innen / inside / interior 550 x 265 x 115 mm	innen / inside / interior 550 x 265 x 130 mm	innen / inside / interior 550 x 265 x 180 mm	innen / inside / interior 550 x 265 x 240mm

Boden ungelocht / Bottom non perforated / Fondo sin perforar



	88.510.10D		88.510.13D		88.510.15D		88.510.20D		88.510.26D
	88.520.10D		88.520.13D		88.520.15D		88.520.20D		88.520.26D
	88.530.10D		88.530.13D		88.530.15D		88.530.20D		88.530.26D
	88.540.10D		88.540.13D		88.540.15D		88.540.20D		88.540.26D
	88.550.10D		88.550.13D		88.550.15D		88.550.20D		88.550.26D
	88.560.10D		88.560.13D		88.560.15D		88.560.20D		88.560.26D

Boden gelocht / Bottom perforated / Fondo perforado



	88.511.10D		88.511.13D		88.511.15D		88.511.20D		88.511.26D
	88.521.10D		88.521.13D		88.521.15D		88.521.20D		88.521.26D
	88.531.10D		88.531.13D		88.531.15D		88.531.20D		88.531.26D
	88.541.10D		88.541.13D		88.541.15D		88.541.20D		88.541.26D
	88.551.10D		88.551.13D		88.551.15D		88.551.20D		88.551.26D
	88.561.10D		88.561.13D		88.561.15D		88.561.20D		88.561.26D

Zubehör

Accessories

Accesorios

Container Deckel

Container lid

Tapa de contenedor

580 x 280 mm

	88.564.01D
	88.564.02D
	88.564.03D
	88.564.04D
	88.564.05D
	88.564.06D



Sicherheitsdeckel

Safety lid

Tapa de seguridad


580 x 280 mm

	88.565.01DK
	88.565.02DK
	88.565.03DK
	88.565.04DK
	88.565.05DK
	88.565.06DK



		Boden ungelocht Bottom non perforated Fondo sin perforar	Boden gelocht Bottom perforated Fondo perforado
außen / outside / exterior	innen / inside / interior		
580 x 280 x 100 mm	550 x 262 x 84 mm	88.566.10D	88.565.10D
580 x 280 x 135 mm	550 x 262 x 117 mm	88.566.13D	88.565.13D
580 x 280 x 150 mm	550 x 262 x 132 mm	88.566.15D	88.565.15D
580 x 280 x 200 mm	550 x 262 x 182 mm	88.566.20D	88.565.20D
580 x 280 x 260 mm	550 x 262 x 244 mm	88.566.26D	88.565.26D





DEWIMED Aluminium 3/4 Container

DEWIMED Aluminium 3/4 Container

DEWIMED Contenedor de Aluminio 3/4



1-1200 **Zyklen/Cycles/Ciclos**

100	135	150
außen / outside / exterior 465 x 280 x 100 mm	außen / outside / exterior 465 x 280 x 135 mm	außen / outside / exterior 465 x 280 x 150 mm
innen / inside / interior 425 x 260 x 84 mm	innen / inside / interior 425 x 260 x 115 mm	innen / inside / interior 425 x 260 x 130 mm

Boden ungelocht / Bottom non perforated / Fondo sin perforar



	88.410.10D		88.410.13D		88.410.15D
	88.420.10D		88.420.13D		88.420.15D
	88.430.10D		88.430.13D		88.430.15D
	88.440.10D		88.440.13D		88.440.15D
	88.450.10D		88.450.13D		88.450.15D
	88.460.10D		88.460.13D		88.460.15D

Boden gelocht / Bottom perforated / Fondo perforado



	88.411.10D		88.411.13D		88.411.15D
	88.421.10D		88.421.13D		88.421.15D
	88.431.10D		88.431.13D		88.431.15D
	88.441.10D		88.441.13D		88.441.15D
	88.451.10D		88.451.13D		88.451.15D
	88.461.10D		88.461.13D		88.461.15D


Zubehör
Accessories
Accesorios

Container Deckel

Container Lid

Tapa de Contenedor

465 x 280 mm

	88.464.01D
	88.464.02D
	88.464.03D
	88.464.04D
	88.464.05D
	88.464.06D




Sicherheitsdeckel

Safety Lid

Tapa de seguridad

465 x 280 mm

	88.464.01DK
	88.464.02DK
	88.464.03DK
	88.464.04DK
	88.464.05DK
	88.464.06DK




		Boden ungelocht Bottom non perforated Fondo sin perforar	Boden gelocht Bottom perforated Fondo perforado
außen / outside / exterior	innen / inside / interior		
465 x 280 x 100 mm	425 x 262 x 85 mm	88.465.10D	88.466.10D
465 x 280 x 135 mm	425 x 262 x 121 mm	88.465.13D	88.466.13D
465 x 280 x 150 mm	425 x 262 x 135 mm	88.465.15D	88.466.15D



DEWIMED Aluminium 1/2 Container

DEWIMED Aluminium 1/2 Container

DEWIMED Contenedor de Aluminio 1/2



1-1200 **Zyklen**/Cycles/Ciclos

100	135	150	200	260
außen / outside / exterior 285 x 280 x 100 mm	außen / outside / exterior 285 x 280 x 135 mm	außen / outside / exterior 285 x 280 x 150 mm	außen / outside / exterior 285 x 280 x 200 mm	außen / outside / exterior 285 x 280 x 260 mm
innen / inside / interior 262 x 260 x 84mm	innen / inside / interior 262 x 260 x 116 mm	innen / inside / interior 262 x 260 x 132 mm	innen / inside / interior 262 x 260 x 182 mm	innen / inside / interior 262 x 260 x 238 mm

Boden ungelocht / Bottom non perforated / Fondo sin perforar



	88.610.10D		88.610.13D		88.610.15D		88.610.20D		88.610.26D
	88.620.10D		88.620.13D		88.620.15D		88.620.20D		88.620.26D
	88.630.10D		88.630.13D		88.630.15D		88.630.20D		88.630.26D
	88.640.10D		88.640.13D		88.640.15D		88.640.20D		88.640.26D
	88.650.10D		88.650.13D		88.650.15D		88.650.20D		88.650.26D
	88.660.10D		88.660.13D		88.660.15D		88.660.20D		88.660.26D

Boden gelocht / Bottom perforated / Fondo perforado



	88.611.10D		88.611.13D		88.611.15D		88.611.20D		88.611.26D
	88.621.10D		88.621.13D		88.621.15D		88.621.20D		88.621.26D
	88.631.10D		88.631.13D		88.631.15D		88.631.20D		88.631.26D
	88.641.10D		88.641.13D		88.641.15D		88.641.20D		88.641.26D
	88.651.10D		88.651.13D		88.651.15D		88.651.20D		88.651.26D
	88.661.10D		88.661.13D		88.661.15D		88.661.20D		88.661.26D

Zubehör

Accessories

Accesorios

Container Deckel

Container Lid

Tapa de Contenedor

285 x 280 mm

	88.655.01D
	88.655.02D
	88.655.03D
	88.655.04D
	88.655.05D
	88.655.06D



Sicherheitsdeckel

Safety Lid

Tapa de Seguridad


285 x 280 mm

	88.656.01DK
	88.656.02DK
	88.656.03DK
	88.656.04DK
	88.656.05DK
	88.656.06DK



		Boden ungelocht Bottom non perforated Fondo sin perforar	Boden gelocht Bottom perforated Fondo perforado
außen / outside / exterior	innen / inside / interior		
285 x 280 x 100 mm	262 x 262 x 84 mm	88.656.10D	88.657.10D
285 x 280 x 135 mm	262 x 262 x 117 mm	88.656.13D	88.657.13D
285 x 280 x 150 mm	262 x 262 x 132 mm	88.656.15D	88.657.15D
285 x 280 x 200 mm	262 x 262 x 182 mm	88.656.20D	88.657.20D
285 x 280 x 260 mm	262 x 262 x 244 mm	88.656.26D	88.657.26D





DEWIMED Aluminium Flachcontainer

DEWIMED Aluminium Flat Container

DEWIMED Contenedor de Aluminio Plano



1-1200 **Zyklen**/Cycles/Ciclos

55	85
außen / outside / exterior 285 x 280 x 55 mm	außen / outside / exterior 285 x 280 x 85 mm
innen / inside / interior 270 x 270 x 33 mm	innen / inside / interior 270 x 270 x 64 mm

Boden ungelocht / Bottom non perforated / Fondo sin perforar



	88.610.05L
	88.620.05L
	88.630.05L
	88.640.05L
	88.650.05L
	88.660.05L

Boden gelocht / Bottom perforated / Fondo perforado



	88.610.08L		88.611.08L
	88.620.08L		88.621.08L
	88.630.08L		88.631.08L
	88.640.08L		88.641.08L
	88.650.08L		88.651.08L
	88.660.08L		88.661.08L

Container Deckel

Container Lid

Tapa de Contenedor


465 x 280 mm

	88.662.01L
	88.662.02L
	88.662.03L
	88.662.04L
	88.662.05L
	88.662.06L



		Boden ungelocht Bottom non perforated Fondo sin perforar	Boden gelocht Bottom perforated Fondo perforado
außen / outside / exterior	innen / inside / interior		
285 x 280 x 55 mm	270 x 270 x 33 mm	88.662.50L	
285 x 280 x 85 mm	270 x 270 x 64 mm	88.662.60L	88.662.65L





DEWIMED Aluminium Mini Container

DEWIMED Aluminium Mini Container

DEWIMED Mini Contenedor de Aluminio



1-1200 **Zyklen/Cycles/Ciclos**

40	70	100
außen / outside / exterior 300 x 140 x 40 mm	außen / outside / exterior 300 x 140 x 70 mm	außen / outside / exterior 300 x 140 x 100 mm
innen / inside / interior 280 x 130 x 27 mm	innen / inside / interior 280 x 130 x 55 mm	innen / inside / interior 280 x 130 x 84 mm

Boden ungelocht / Bottom non perforated / Fondo sin perforar



	88.660.04L		88.660.07L		88.660.10L
	88.665.04L		88.665.07L		88.665.10L
	88.670.04L		88.670.07L		88.670.10L
	88.675.04L		88.675.07L		88.675.10L
	88.680.04L		88.680.07L		88.680.10L
	88.685.04L		88.685.07L		88.685.10L

Boden gelocht / Bottom perforated / Fondo perforado



	88.661.04L		88.661.07L		88.661.10L
	88.666.04L		88.666.07L		88.666.10L
	88.671.04L		88.671.07L		88.671.10L
	88.676.04L		88.676.07L		88.676.10L
	88.681.04L		88.681.07L		88.681.10L
	88.686.04L		88.686.07L		88.686.10L

Zubehör

Accessories

Accesorios

Container Deckel

Container Lid

Tapa de Contenedor


465 x 280 mm

	88.687.01L
	88.687.02L
	88.687.03L
	88.687.04L
	88.687.05L
	88.687.06L



		Boden ungelocht Bottom non perforated Fondo sin perforar	Boden gelocht Bottom perforated Fondo perforado
außen / outside / exterior	innen / inside / interior		
300 x 140 x 40 mm	280 x 130 x 27 mm	88.687.50L	88.687.55L
300 x 140 x 70 mm	280 x 130 x 55 mm	88.687.60L	88.687.65L
300 x 140 x 100 mm	280 x 130 x 84 mm	88.687.70L	88.687.75L






Siebkörbe aus Edelstahl MAX 1/1 Wire Baskets and Trays MAX 1/1 Bandejas de Acero MAX 1/1


Komplettsystem mit Deckel
Complete System with Lid
Sistema Completo con Tapa

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.567.35D	545x255x35 mm
88.567.55D	545x255x55 mm
88.567.75D	545x255x75 mm
88.567.15D	545x255x105 mm
88.568.35D	485x255x35 mm
88.568.55D	485x255x55 mm
88.568.75D	485x255x75 mm
88.568.15D	485x255x105 mm




Deckel für Siebkorb
Lid for Stainless Steel Basket
Tapas para Cesto de Esterilización

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.567.00D	545 x 255 mm
88.568.00D	485 x 255 mm



Siebkorb ohne Deckel
Stainless Steel Basket without Lid
Cesto de Esterilización sin Tapa

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.567.30D	545x255x30 mm
88.567.50D	545x255x50 mm
88.567.70D	545x255x70 mm
88.567.10D	545x255x100 mm
88.568.30D	485x255x30 mm
88.568.50D	485x255x50 mm
88.568.70D	485x255x70 mm
88.568.10D	485x255x100 mm



Siebkörbe aus Edelstahl MAX 3/4 Wire Baskets and Trays MAX 3/4 Bandejas de Acero Inoxidable MAX 3/4


Komplettsystem mit Deckel
Complete System with Lid
Sistema Completo con Tapas

3/4	
Art. No.	Größe - Dimensions - Dimensiones
88.468.35D	405 x 255 x 35mm
88.468.55D	405 x 255 x 55mm
88.468.75D	405 x 255 x 75mm
88.468.15D	405 x 255 x 105mm



Deckel für Siebkorb
Lid for Stainless Steel Basket
Tapa para Cesto de Esterilización


3/4	
Art. No.	Größe - Dimensions - Dimensiones
88.468.00D	405 x 255 mm



Siebkorb ohne Deckel
Stainless Steel Basket without Lid
Cesto de Esterilización sin Tapas

3/4	
Art. No.	Größe - Dimensions - Dimensiones
88.468.30D	405 x 255 x 30mm
88.468.50D	405 x 255 x 50mm
88.468.70D	405 x 255 x 70mm
88.468.10D	405 x 255 x 100mm






Siebkörbe aus Edelstahl MAX 1/2 Wire baskets and Trays MAX 1/2 Bandejas de Acero MAX 1/2


Komplettsystem mit Deckel
Complete System with Lid
Sistema Completo con Tapa

1/2	
Art. No.	Größe - Dimensions - Dimensiones
88.658.35D	255 x 245 x 35 mm
88.658.55D	255 x 245 x 55 mm
88.658.70D	255 x 245 x 70 mm
88.658.15D	255 x 245 x 105 mm




Deckel für Siebkorb
Lid for Stainless Steel Basket
Tapa para Cesto de Esterilización

1/2	
Art. No.	Größe - Dimensions - Dimensiones
88.658.00D	255 x 245 mm



Siebkorb ohne Deckel
Stainless Steel Basket without Lid
Cesto de Esterilización sin Tapas

1/2	
Art. No.	Größe - Dimensions - Dimensiones
88.658.30D	255 x 245 x 30 mm
88.658.50D	255 x 245 x 50 mm
88.658.65D	255 x 245 x 65 mm
88.658.10D	255 x 245 x 100 mm



Siebkörbe aus Edelstahl

Wire Baskets and Trays

Bandeja de Acero

Siebkorb ohne Deckel
Stainless Steel Basket without Lid
Cesto de Esterilización sin Tapas

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.806.03D	484x252x25mm
88.806.05D	484x252x45mm
88.806.07D	484x252x65mm
88.806.10D	484x252x95mm
88.807.03D	539x252x25mm
88.807.05D	539x252x45mm
88.807.07D	539x252x65mm
88.807.10D	539x252x95mm

3/4	
Art. No.	Größe - Dimensions - Dimensiones
88.805.03D	405x252x25mm
88.805.05D	405x252x45mm
88.805.07D	405x252x65mm
88.805.10D	405x252x95mm


1/2	
Art. No.	Größe - Dimensions - Dimensiones
88.804.03D	242x252x25mm
88.804.05D	242x252x45mm
88.804.07D	242x252x65mm
88.804.10D	242x252x95mm

Deckel für Siebkorb
Lid for Stainless Steel Basket
Tapas para Cesto de Esterilización

1/1, 3/4, 1/2	
Art. No.	Größe - Dimensions - Dimensiones
88.812.01D	246x256x15mm
88.812.02D	409x256x15mm
88.812.03D	488x256x15mm
88.812.04D	543x256x15mm

Art. No.	Größe - Dimensions - Dimensiones
88.895.03D	Füße schwarz, aus Kunststoff/ Stands, black, made of plastic/ Soporte, negro, hecho de plástico





Drahtsiebkörbe aus Edelstahl **Wire baskets made from Stainless Steel** **Bandeja de Acero Inoxidable**


Drahtsiebkorb Edelstahl
Stainless Steel wire Basket
Bandeja de Acero Inoxidable

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.591.03	540x255x30 mm
88.591.05	540x255x50 mm
88.591.07	540x255x70 mm
88.591.10	540x255x100 mm
88.590.03	485x255x30 mm
88.590.05	485x255x50 mm
88.590.07	485x255x70 mm
88.590.10	485x255x100 mm
3/4	
88.467.03	405x255x30 mm
88.467.05	405x255x50 mm
88.467.07	405x255x70 mm
88.467.10	405x255x100 mm
1/2	
88.690.03	255x245x30 mm
88.690.05	255x245x50 mm
88.690.07	255x245x70 mm
88.690.10	255x245x100 mm




Zwischenboden
Wire Base
Fondo Intermedio


1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.598.30	545x255x30 mm
3/4	
88.468.03	405x255x30 mm
1/2	
88.699.30	258x258x30 mm



Spannrahmen
Drape Retainer
Fijador de Pano

1/1	
Art. No.	Größe - Dimensions - Dimensiones
88.595.00	545x255 mm
3/4	
88.468.05	405x255 mm
1/2	
88.695.00	255x255 mm






Siebkörbe aus Edelstahl für Mini Container

Baskets made from Stainless Steel for Mini Container

Bandeja de Acero Inoxidable para Mini Contenedor


Drahtsiebkorb Edelstahl
Stainless Steel wire Basket
Bandeja de Acero Inoxidable

Mini	
Art. No.	Größe - Dimensions - Dimensiones
88.784.45D	265 x 120 x 45 mm




Siebkorb aus Edelstahl
Basket made from Stainless Steel
Bandeja de Acero Inoxidable

Mini	
Art. No.	Größe - Dimensions - Dimensiones
88.785.22D	268 x 124 x 22mm



Siebkorb aus Edelstahl mit / ohne Deckel
Basket made from Stainless Steel with / without Lid
Bandeja de Acero Inoxidable con / sin Tapa


Mini	
Art. No.	Größe - Dimensions - Dimensiones
Siebkorb mit Deckel Stainless Steel Basket with Lid Cesto de Esterilización con Tapa	
88.785.45D	271 x 117 x 45 mm
Siebkorb ohne Deckel Stainless Steel Basket without Lid Cesto de Esterilización sin Tapa	
88.785.40D	271 x 117 x 40 mm
Deckel für Siebkorb Lid for Stainless Steel Basket Tapas para Cesto de Esterilización	
88.785.00D	271 x 117 mm



Verschiedene Siebkörbe Varieties of Baskets Variedad de Bandejas


Ophthalmologie-Siebkorb mit Silikonhalterung Ophthalmic Basket with Silicon Holder Bandeja de Oftalmología con Silicona

Art. No.	Größe - Dimensions - Dimensiones
88.795.24D	240 x 250 x 38 mm
88.795.48D	480 x 250 x 38 mm



Micro- Siebkorb mit Silikonhalterung Micro Basket with Silicon Holder Micro Bandeja con Silicona


Art. No.	Größe - Dimensions - Dimensiones
88.795.34D	340 x 250 x 40 mm
88.795.44D	440 x 250 x 40 mm



Siebkorb mit feinem Drahtgewebe Basket with fine Wire Mesh Bandeja con Alambre Fino

Art. No.	Größe - Dimensions - Dimensiones
88.796.12D	120 x 155 x 40 mm
88.796.23D	235 x 155 x 40 mm
88.796.38D	380 x 270 x 40 mm
88.796.17D	275 x 175 x 35 mm
88.796.20D	275 x 200 x 35 mm





Kleiner Siebkorb

Small Basket

Caja Pequeña

Art. No.	Größe - Dimensions - Dimensiones
88.797.90D	90 x 60 x 25 mm
88.797.10D	105 x 70 x 25 mm
88.797.12D	120 x 80 x 25 mm



Endoskopie Siebkorb für 1 Endoskop

Endoscopy Basket for 1 Endoscope

Caja de Endoscopia para 1 Endoscopio

Art. No.	Größe - Dimensions - Dimensiones
88.798.44D	272 x 83 x 44 mm



Endoskopie Siebkorb für 2 Endoskope

Endoscopy Basket for 2 Endoscopes

Caja de Endoscopia para 2 Endoscopios

Art. No.	Größe - Dimensions - Dimensiones
88.799.44D	272 x 103 x 44 mm



Endoskopie Siebkorb für 1 Endoskop

Endoscopy Basket for 1 Endoscope

Caja de Endoscopia para 1 Endoscopio

Art. No.	Größe - Dimensions - Dimensiones
88.800.44D	452 x 83 x 44 mm


Endoskopie Siebkorb für 2 Endoskope

Endoscopy Basket for 2 Endoscopes

Caja de Endoscopia para 2 Endoscopios


Art. No.	Größe - Dimensions - Dimensiones
88.801.44D	452 x 103 x 44 mm



Endoskopie Rack

Endoscopy Rack

Rack de Endoscopia


Art. No.	Größe - Dimensions - Dimensiones
88.800.00D	485 x 235 x 160 mm






Zubehör Accessories Accesorios

Rundfilter für 1/1, 3/4 und 1/2 Container (ø 190 mm)
Round filter for 1/1, 3/4 and 1/2 Container (ø 190 mm)
Filtro redondo para 1/1, 3/4 y 1/2 Contenedor (ø 190 mm)



Einweg Papier Filter Disposable Paper Filter Papel Disponible	Textil Filter (50 Zyklen) Textile Filter (50 cycle) Filtro Textil (50 ciclos)	PTFE Filter (1200 Zyklen) PTFE Filter (1200 cycle) PTFE Filtro (1200 ciclos)
Art. No. 88.700.00	Art. No. 88.700.01	Art. No. 88.701.00
1000 pcs.	50 pcs.	1 pc.

Filter für Mini und Dental Container
Filter for Mini and Dental Container
Filtro para Contenedor Mini y Dental




Einweg Papier Filter Disposable Paper Filter Papel Disponible	Textil Filter (50 Zyklen) Textile Filter (50 cycle) Filtro Textil (50 ciclos)	PTFE Filter (1200 Zyklen) PTFE Filter (1200 cycle) PTFE Filtro (1200 ciclos)
Art. No. 88.700.04	Art. No. 88.700.05	Art. No. 88.701.01
1000 pcs.	50 pcs.	1 pc.

Filterhalter
Filter holding device
Filtro que Sostiene Dispositivo



für 1/1, 3/4, 1/2 und Flachcontainer for 1/1, 3/4, 1/2 and Flat Container para 1/1, 3/4, 1/2 y Contenedor Plano	für Mini und Dental Container for Mini and Dental Container para Contenedor Mini y Dental
Art. No.	Art. No.
88.700.10	88.700.20

Beschriftungsschilder mit Indikator
Paper Labels with Indicator
Etiqueta de Papel con Indicador




Groß (67 x 40 mm) Large (67 x 40 mm) Grande (67 x 40 mm)	Klein (60 x 18 mm) Small (60 x 18 mm) Chico (60 x 18 mm)
88.710.02	88.710.04
1000 pcs.	1000 pcs.

Sicherheitsplombe aus Plastik
Plastic Security Seal
Precinto de Plástico



Art. No. 88.710.10




Identifikationsschilder (Max. 12 Zeichen) Identification Tags (Max. 12 Letters) Placa de Identificaciones (Max. 12 letras)

Schild ohne Beschriftung Tag without Marking Placa sin Etiquetado	
	88.710.20
	88.710.22
	88.710.24
	88.710.26
	88.710.28
	88.710.29

Schild mit Beschriftung Tag with Marking Placa con Etiquetado	
	88.710.20-TXT
	88.710.22-TXT
	88.710.24-TXT
	88.710.26-TXT
	88.710.28-TXT
	88.710.29-TXT

Schild mit Beschriftung incl. Haken Tag with Marking and Hook Placa con Etiquetado y Gancho	
	88.710.20-TXT
	88.710.22-TXT
	88.710.24-TXT
	88.710.26-TXT
	88.710.28-TXT
	88.710.29-TXT



Unterteilung und Sortierstifte Partition Sheets and Holding Pins Separadores y Pin para Instrumentos

Sortierstift für Instrumente Holding Pin for Instruments Pin para Instrumentos	Sortierstift mit Schlitz Holding Pin for Partition Sheets Pin para Instrumentos Ranurado	Trennbleche in mm Sheets in mm Separadores en mm	
Art. 88.720.00	88.720.01	50 x 20	88.720.05
		130 x 20	88.720.13
		225 x 20	88.720.22
		460 x 20	88.720.46


Halteklammer
Fixation Clamp
Clamp de Fijación

Einfach / Single / Sencillo	Art. No. 88.321.01
Doppelt / Double / Doble	Art. No. 88.321.02




Sterilisationstuch
Sterilisation Drape
Paño de Esterilización


Abmessungen / Size / Tamaño	Grün / Green / Verde
40 x 60 cm	88.260.04
100 x 100 cm	88.260.10
100 x 140 cm	88.260.13
100 x 160 cm	88.260.16



Silikonmatte, gelocht
Silicon Mat, Perforated
Alfombra de Silicona, Perforada

Abmessungen / Size / Tamaño	Blau / Blue / Azul
270 x 125 mm	88.722.27
245 x 245 mm	88.722.24
220 x 230 mm	88.722.22
380 x 230 mm	88.722.38
440 x 230 mm	88.722.44
520 x 230 mm	88.722.52
340 x 100 mm	88.722.34







DEWIMED Medizintechnik GmbH

Unter Hasslen 14

D-78532 Tuttlingen

Germany

+49-07462-923 930

+49-07462-923 9333

info@dewimed.de

Cat No REF. 00.088.10 Ver. 02/2014

www.dewimed.de

Graphic Design: Hichem Boukerch