




NO SCALPEL VASECTOMY INSTRUMENTS

SET DE VASECTOMIA
NO INCISIVO




1 Fixation of the vas under the scrotal skin using the extracutaneous vas fixation clamp DEWIMED no-scalpel vasectomy fixator ring clamp




1 Fijación del conducto bajo la piel del escroto utilizando el clamp DEWIMED de fijación externo para vasectomía sin bisturí

2 Puncture of the scrotal skin using the DEWIMED No-Scalpel Vasectomy forceps with sharp pointed smooth curved jaw




2 Punción de la piel de escroto usando la pinza de vasectomía DEWIMED sin bisturí aguda cortante y curva

3 Dissecting blood vessels and sheath of vas prior to division and ligation or cauterization of thlumen



3 Dirección de los vasos sanguíneos y del conducto antes de la división y ligadura o cauterización del lumen



Set

No. de Catálogo	Piezas	Descripción
12.416.14	2	Pinza para vasectomía de anillo de ligación medida: 14 cm.
12.417.14	1	Pinza para vasectomía de ligación con punta aguda medida: 14 cm.
24.148.13	1	BROWN porta agujas medida: 13 cm.
12.221.14	2	HALSTED mosquito pinza medida: 14 cm.
08.281.15	1	METZENBAUM tijeras curvas medida: 15 cm.
10.180.15	1	ADSON pinza 1 x 2 dientes medida: 15 cm.
18.245.01	1	Gancho para vasectomía medida: 15 cm.
88.676.07	1	Contenedor especial (tapa y fondo perforados) medida: 300 x 140 x 70 mm



Set

No. of Catalog	Pieces	Description
12.416.14	2	Vasectomy forceps with blunt ring tips measure: 14 cm.
12.417.14	1	Vasectomy forceps with sharp pointed smooth curved jam measure: 14 cm.
24.148.13	1	BROWN needle holders measure: 13 cm.
12.221.14	2	HALSTED mosquito forceps measure: 14 cm.
08.281.15	1	METZENBAUM curved scissors measure: 15 cm.
10.180.15	1	ADSON forceps 1 x 2 teeth measure: 15 cm.
18.245.01	1	Vasectomy hook measure: 15 cm.
88.676.07	1	Special container (lid and bottom perforated) measure: 300 x 140 x 70 mm



DEWIMED Medizintechnik GmbH
Unter Hasslen 14
D-78532 Tuttlingen
Germany

+49-07462-923 930
+49-07462-923 9333
@ info@dewimed.de

Technische Änderungen vorbehalten - Kein Teil des Kataloges darf in irgend einer Form ohne schriftliche Genehmigung unseres Hauses reproduziert, vervielfältigt oder verbreitet werden.

We reserve the right to make alterations - Any part of this catalogue may not be reproduced, copied or distributed in any way without our written permission.

Cambios técnicos reservados - Se prohíbe copiar, reproducir o distribuir cualquier parte de este catálogo sin nuestro permiso por escrito.

Printed in Mexico DEWIMED Medizintechnik GmbH D-78532 Tuttlingen.